

DIANE

A BOW MEDICAL brand

BM-DEV-EN-017 - Release note

Version 4.14.3.2

OS Supported: Win10, Win11, WinServer2016, WinServer2019, WinServer2022. Vidal API Mini: 2020.5.3
Vidal API Rec.: 2024.11.00 BCB API Mini.: 20180702

CORRECTIONS

Desktop client

DKC-63473

Fixed DiaUpdate error on new installation

DKC-63478 - (Associated with the medical device)

Correction of a problem updating the care plan when linking an occurrence to a PSE on a free-flow prescription occurrence.

DKC-63509 - (Associated with the medical device)

Correction of a problem with the resumption of a PSE link following a relay.